

## Azura ExtremeTaq DNA Polymerase

Catalog No.	Pack Size and Concentration	Components and Volume	
AZ-1802	250 units, 5u/μl	ExtremeTaq Polymerase - 1 x 50 μl	5x ExtremeTaq Buffer – 2 x 1ml
AZ-1810	1000 units, 5u/μl	ExtremeTaq Polymerase - 4 x 50 μl	5x ExtremeTaq Buffer – 8 x 1ml

### Description

Azura ExtremeTaq DNA Polymerase is a high-performance DNA polymerase complex comprised of Azura HS Taq and a proof-reading polymerase in a unique balance which is ideally suited to problematic GC-rich and AT-rich genomic DNA amplification. Azura ExtremeTaq generates extremely high DNA yields and 3'-A overhangs which are ideal for direct integration into TA cloning vectors. In addition, the enzyme delivers exceptional PCR sensitivity in low copy assays with 4-fold higher fidelity than wild-type *Taq* polymerase. The DNA Polymerase is provided in a simple two-component format for ease of use, minimal handling and greater reproducibility. Azura ExtremeTaq DNA Polymerase is supplied with a complete 5x reaction buffer inclusive of GC enhancers, optimal levels of dNTP and MgCl<sub>2</sub>, and a high ionic strength for broad compatibility with amplicon size and complexity.

- New-generation 5x PCR buffer formulation including optimal levels of dNTPs, MgCl<sub>2</sub> and PCR enhancers for maximum PCR efficiency and reaction speed.
- Robust PCR performance across a wide range of DNA templates including low-copy assays and problematic, long templates.
- High-yields with amplicons up to 10 Kb with standard or fast cycling.

### Storage

Azura ExtremeTaq DNA Polymerase is shipped on blue or dry ice and should be stored at –20°C upon receipt. Excessive freeze/thawing should be avoided. When stored as specified, Azura ExtremeTaq DNA Polymerase is stable for 12 months from date of receipt. The Kit may also be stored at 4°C for 1 month.

### Important Guidelines

**Azura 5x ExtremeTaq Buffer:** The 5x reaction buffer contains proprietary PCR enhancers, optimal levels of dNTPs (5mM) and 15mM MgCl<sub>2</sub>. The buffer has been designed to deliver maximum efficiency, sensitivity and success with difficult amplicons. We do not suggest the use of additional PCR enhancers or dNTPs.

**Template:** For complex genomic DNA, we suggest the use of 5ng - 500ng per reaction; For cDNA or plasmid DNA, please use < 100ng per reaction.

**Primers:** Primers should have a predicted melting temperature of around 60°C, using default Primer 3 settings (<http://frodo.wi.mit.edu/primer3/>). The final primer concentration in the reaction should be between 0.2μM and 0.6μM.

**Annealing:** We recommend performing a temperature gradient to determine the optimal annealing temperature. Alternatively, 57°C can be used as a starting point. For optimization, increase in 2°C increments. For example, if non-specific products are present or smearing is visible, a higher annealing temperature is required.

**Extension:** Optimal extension is achieved at 72°C. The optimal extension time is dependent on amplicon length and complexity.

15 seconds per kilobase(kb) is recommended for amplification from eukaryotic genomic DNA or cDNA amplicons less than 5kb. For amplicons which are greater than 5kb, we suggest 40 seconds to 60 seconds per kb.

## Reaction setup

1. Prepare a PCR master mix based on following table:

Component	50µl Reaction	Final Concentration/Notes
Azura 5x ExtremeTaq Buffer	10 µl	1x
Forward Primer (10µM)	2.0 µl	400 nM
Reverse Primer (10µM)	2.0 µl	400 nM
Template DNA	<100ng cDNA, <500ng genomic DNA	variable
Azura ExtremeTaq DNA Polymerase (5u/µl)	0.25 µl - 1 µl	variable
PCR-grade water	Up to 50 µl final volume	

\* For alternative total reaction volumes (eg. 20 µl), scale all components proportionally to maintain final concentrations.

2. PCR cycling:

Cycles	Temperature & Time	Notes
1	95°C, 1 minute to 2 minutes	Initial Denaturation, enzyme activation
25 - 40	95°C, 15 seconds 55°C to 67°C, 15 seconds 72°C, 15 seconds per Kb*	Denaturation Annealing (determined by user) Extension* (see notes above)

## Quality Control

Azura ExtremeTaq DNA Polymerase is tested extensively for robust activity, processivity, efficiency, heat activation, sensitivity, absence of nuclease contamination and absence of nucleic acid contamination. Azura ExtremeTaq DNA Polymerase is manufactured under a comprehensive quality management system, following ISO 9001:2008 standards.

## Limitations of Use

This product is intended for research purposes only and is not intended for any animal or human therapeutic use.

## Technical Support

For Trouble-shooting and Technical Guidance, please contact us at [tech@azuragenomics.com](mailto:tech@azuragenomics.com) and provide PCR reaction conditions, cycling parameters, amplicon size, and screen grabs (gel images) if possible.

Azura Genomics guarantees the performance of all products in the manner described in our product documentation. The purchaser must determine the suitability of the product for its particular use. Azura Genomics shall not be liable for any direct, indirect, consequential or incidental damages arising from the use, the result of use or the liability to use this product.