

## Proteinase K

| Catalog No. | Pack Size | Components                |
|-------------|-----------|---------------------------|
| AZ-3101     | 100 mg    | Proteinase K - 1 x 100 mg |
| AZ-3105     | 500 mg    | Proteinase K - 5 x 100 mg |
| AZ-3110     | 1000 mg   | Proteinase K- 10 x 100 mg |

### Description

Proteinase K is provided as a lyophilized powder form and is of molecular biology grade. It is a broad-spectrum serine protease with very high specific activity. Proteinase K is a recombinant enzyme expressed in *Pichia pastoris*, and undergoes extensive purification to yield the highest quality product. It is active under a wide range of reaction conditions, including elevated temperatures and in the presence of SDS. As a result, Proteinase K is widely using for the digestion of proteins, including DNases and RNases, during nucleic acid preparations without compromising the integrity of isolated DNA or RNA. Proteinase K is free of exonucleases, endonucleases, and ribonucleases.

- Recombinant broad-spectrum non-specific protease derived from *Tritirachium album* and over-expressed in *Pichia pastoris*.
- Molecular Biology Grade
- High activity and exceptional purity.
- Active at high temperatures (up to 56 °C) and denaturing conditions (in the presence of urea and/or SDS), making it ideal for digesting proteins in variety of applications.
- Stable over a wide pH range: 4.0–12.5 (optimum pH 7.5–8.0).
- Low residual DNA content ( $\leq 10$  pg/mg).

### Storage

Storage of Proteinase K at -20°C is recommended.

### Recommendations for Use

Proteinase K is soluble in water, PBS and Tris. We recommend dissolving our lyophilized Proteinase K powder into a 20 mg/mL solution of 50 mM Tris-HCl of pH = 7.8 and 3 mM CaCl<sub>2</sub> for immediate use, which gives an activity of  $\geq 800$  U/ml. However, if you intend to prepare Proteinase K solution for long-term storage -20°C, dissolve the powder in 50 mM Tris-HCl of pH = 7.8, 3 mM CaCl<sub>2</sub> and 50% glycerol.

### Important Guidelines

**Protein concentration:** Protein concentration is determined by measuring absorbance at 280 nm.

**Exonuclease activity:** Free of detectable exonucleases activity as judged by gel electrophoresis following incubation of 1 µg of HindIII-digested λ DNA with 50 µg of Proteinase K for 16 h at 37°C.

**Endonuclease activity:** Free of detectable endonucleases activity as judged by gel electrophoresis following incubation of 1 µg pUC19 DNA with 40 µg of Proteinase K for 16 h at 37°C.

**RNase activity:** Free of detectable RNase activities as judged by gel electrophoresis following incubation of 2 µg rRNA from E. coli with 20 µg of Proteinase K for 4 h at 37°C.

**DNA content:** DNA content is  $\leq 10$  pg/mg, which is determined by qPCR.

## Quality Control

Our Proteinase K activity measurements demonstrate very low batch-to-batch variability. Such a high reproducibility enables stable working conditions, and therefore repeatable and reliable experiment results.

|                            | <b>Proteinase K- Lyophilized Powder, MB Grade</b>   |
|----------------------------|---|
| <b>Source</b>              | <i>Parengyodontium album (Tritirachium album)</i>   |
| <b>Host</b>                | <i>Komagataella phaffii (Pichia pastoris)</i>   |
| <b>Solubility in Water</b> | ≥ 20 mg/ml  |
| <b>Activity</b>            | ≥ 30 U/mg lyophilizate  |
| <b>Specific Activity</b>   | ≥ 40 U/mg protein   |
| <b>Unit Definition</b>     | One unit of Proteinase K hydrolyzes urea-denatured hemoglobin producing color equivalent of 1 µmol tyrosine per 1 min at 37°C and pH 7.5 (Folin & Ciocalteu's method), 1 U = 1 mAnsonU. |
| <b>Protein Content</b>     | ≥ 70%   |
| <b>DNA Contamination</b>   | ≤ 10 pg/mg  |
| <b>Storage Conditions</b>  | -20°C   |
| <b>Shelf Life</b>          | 24 months.  |
|                            |   |

## Limitations of Use

This product is intended for research purposes and is not intended for any animal or human therapeutic use.

## Technical Support

For trouble-shooting and technical guidance, please contact us at [tech@azuragenomics.com](mailto:tech@azuragenomics.com) and provide reaction parameters.